
- 1) Detector section (4 detectors)
Analyzes samples. Each detector is equipped with a sensor to detect reagent addition and automatically trigger mixing and optical clot detection.
- 2) Pipette guide
Used to hold the pipette chip when reagent is added and to cover the detector. The pipette guide can be turned in two directions.
- 3) Detector LED
Indicates the timing to add reagent and the instrument's status.
- 4) Analysis start switch
Press this button to analyze sample.
- 5) Incubator
Warms the reagent in the reagent cup.
- 6) Reaction tube holder
Used to set reaction tubes for ease of sample pipetting.
- 7) Built-in printer
Prints analysis results and analyzer settings.
- 8) Input key
For operation and settings in the CA-50 menu.
- 9) LCD screen
Displays analysis conditions, analysis results, and error messages.

<RESULT DISPLAY>

When replication mode is selected, the mean of two results obtained from two individual detector channel is displayed as a test result. Original results can be seen in the "Detailed Result Display".

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